

# Technical Data


Light and compact  
**MICROWAVE  
PRESSURE DIGESTION  
SYSTEM**



- Rapid
- Safe
- Accurate
- Compact



**speed**wave **MWS-2**



# BERGHOF **speedwave** MWS-2

## The advantages of using the **speedwave** MWS-2

- Light, compact stainless-steel oven with a capacity of 27 liters
- Contact-free IR temperature measurement in all pressure vessels
- Very simple, tool-less handling of vessels
- Closeable rotors, ensuring a high degree of safety and long service life
- Integrated suction feature
- Latest Inverter microwave technology
- Integrated control system

## Technical data of the **speedwave** MWS-2

Power supply	230 V / 50 Hz / 1 350 W
Microwave output	1 000 W
Frequency	2 450 MHz
Outer dimensions (WxDxH)	520 x 460 x 330 mm
Interior dimensions (WxDxH)	350 x 340 x 215 mm
Weight	approx. 14 kg
Noise level	< 60 dB
Ambient conditions	15 - 35°C / 85% rel. humidity
Interface	serial (RS-232)

## Digestion vessels for the **speedwave** MWS-2

Vessel type	Volume	T-max. (continuous)	Working pressure	Multiple arrangement
DAP-60K	60 ml	220°C	40 bar / 580 psi	10 vessels
DAC-70	70 ml	220°C	100 bar / 1450 psi	5 vessels
DAQ-10	10 ml	220°C	75 bar / 1090 psi	24 vessels

*Special sizes can also be produced to suit customer applications.*



Ideen  
werden  
wahr.

Presented by:

PRODUCTS + INSTRUMENTS

Ein Unternehmen der BERGHOF Gruppe

### BERGHOF

Products + Instruments GmbH

Harretstrasse 1 • 72800 Eningen • Germany

info@berghof-instruments.de

http://www.berghof-instruments.de

Tel.:+49 7121 894-0 • Fax:+49 7121 894-300